

Mouse IgG2a Isotype Control [MOPC 173] {FITC} (STJA0022659)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	0.5 mg/mL
Conjugation	FITC
Formulation	0.01M PBS pH 7.4, 50% glycerol, 0.05% proclin300
Isotype	IgG2a
Storage Instruction	Store at 4°C for 12 months Protect from light Do not freeze

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS
------	--