

Human VHH-hFc Isotype Control [cAb-Lys2] {FITC} (STJA0022691)

GENERAL INFORMATION

| | |
|---------------|--------------------|
| Product Type | Primary antibodies |
| Applications | FC |
| Host / Source | Human |

PRODUCT PROPERTIES

| | |
|---------------------|---|
| Clonality | Monoclonal |
| Concentration | 0.5 mg/mL |
| Conjugation | FITC |
| Formulation | 0.01M PBS pH 7.4, 50% glycerol, 0.05% proclin300 |
| Isotype | VHH |
| Storage Instruction | Store at 4°C for 12 months Protect from light Do not freeze |

ADDITIONAL INFORMATION

| | |
|------|--|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS |
|------|--|