

## Rat IgG2a Isotype Control [SPE-7] {APC} (STJA0022701)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | FC                 |
| Host / Source | Rat                |

### PRODUCT PROPERTIES

|                     |                                                             |
|---------------------|-------------------------------------------------------------|
| Clonality           | Monoclonal                                                  |
| Concentration       | 0.5 mg/mL                                                   |
| Conjugation         | APC                                                         |
| Formulation         | 0.01M PBS pH 7.4, 50% glycerol, 0.05% proclin300            |
| Isotype             | IgG2a                                                       |
| Storage Instruction | Store at 4°C for 12 months Protect from light Do not freeze |

### ADDITIONAL INFORMATION

|      |                                                                                                                    |
|------|--------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS |
|------|--------------------------------------------------------------------------------------------------------------------|