

Dog IgG2 Isotype Control [HyHEL-10] {APC} (STJA0022720)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	FC
Host / Source	Dog

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	0.5 mg/mL
Conjugation	APC
Formulation	0.01M PBS pH 7.4, 50% glycerol, 0.05% proclin300
Isotype	IgG2
Storage Instruction	Store at 4°C for 12 months Protect from light Do not freeze

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS
------	--