

Human VHH-hFc Isotype Control [cAb-Lys2] {PerCP} (STJA0022779)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	0.5 mg/mL
Conjugation	PerCP
Formulation	0.01M PBS pH 7.4, 50% glycerol, 0.05% proclin300
Isotype	VHH
Storage Instruction	Store at 4°C for 12 months Protect from light Do not freeze

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS
------	--