

**Guinea pig IgG Isotype Control [NA] {PE} (STJA0022787)****GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	FC/IF
Host / Source	Guinea pig

PRODUCT PROPERTIES

Clonality	Polyclonal
Concentration	0.5 mg/mL
Conjugation	PE
Formulation	0.01M PBS pH 7.4, 50% glycerol, 0.05% proclin300
Isotype	IgG
Storage Instruction	Store at 4°C for 12 months Protect from light Do not freeze

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS
------	--