

Feline IgG2 Isotype Control [HyHEL-10] {Unconjugated} (STJA0022829)

GENERAL INFORMATION

| | |
|---------------|--------------------|
| Product Type | Primary antibodies |
| Host / Source | Feline |

PRODUCT PROPERTIES

| | |
|---------------------|---|
| Clonality | Monoclonal |
| Concentration | 0.7 mg/mL |
| Conjugation | Unconjugated |
| Purification | Affinity chromatography |
| Formulation | 0.01M PBS, pH 7.4. |
| Isotype | IgG2 |
| Storage Instruction | Suitable for storage at +4°C between 1-2 weeks For longer term store at -20°C for up to 12 months |

ADDITIONAL INFORMATION

| | |
|------|--|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS |
|------|--|