

Feline IgA Isotype Control [HyHEL-10] {Unconjugated} (STJA0022830)

GENERAL INFORMATION

Product Type	Primary antibodies
Host / Source	Feline

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Affinity chromatography
Formulation	0.01M PBS, pH 7.4.
Isotype	IgA
Storage Instruction	Suitable for storage at +4°C between 1-2 weeks For longer term store at -20°C for up to 12 months

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS
------	--