

## Feline IgE Isotype Control [HyHEL-10] {APC} (STJA0022841)

### GENERAL INFORMATION

Product Type	Primary antibodies
Host / Source	Feline

### PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	APC
Purification	Affinity chromatography
Formulation	0.01M PBS pH 7.4, 50% glycerol, 0.05% proclin300
Isotype	IgE
Storage Instruction	Store at 4°C for 12 months Protect from light Do not freeze

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO) MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS
------	--