

## Anti-Recombinant-PDE2A antibody [RM1C87] (STJA0022883)

| GENERAL INFORMATION    |                     |                                                                                                                      |
|------------------------|---------------------|----------------------------------------------------------------------------------------------------------------------|
|                        | Product Type        | Primary antibodies                                                                                                   |
|                        | Applications        | WB                                                                                                                   |
|                        | Host / Source       | Rabbit                                                                                                               |
|                        | Reactivity          | Human/Mouse/Rat                                                                                                      |
| PRODUCT PROPERTIES     |                     |                                                                                                                      |
|                        | Clonality           | Monoclonal                                                                                                           |
|                        | Conjugation         | Unconjugated                                                                                                         |
|                        | Purification        | Protein A/G purified from cell culture supernatant                                                                   |
|                        | Dilution Range      | WB: 1:1000-1:2000                                                                                                    |
|                        | Formulation         | 0.01M PBS                                                                                                            |
|                        | Isotype             | IgG                                                                                                                  |
|                        | Storage Instruction | Suitable for storage at +4°C between 1-2 weeks. For longer term store at -20°C for up to 12 months.                  |
| TARGET INFORMATION     |                     |                                                                                                                      |
|                        | Gene ID             | <a href="#">5138</a>                                                                                                 |
|                        | Gene Symbol         | <a href="#">PDE2A</a>                                                                                                |
|                        | UniProt ID          | <a href="#">PDE2A_HUMAN</a>                                                                                          |
| ADDITIONAL INFORMATION |                     |                                                                                                                      |
|                        | Note                | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |