

Anti-Recombinant-Bim/BCL2L11 antibody [RM1E10] (STJA0023004)

GENERAL INFORMATION

| | |
|---------------|--------------------|
| Product Type | Primary antibodies |
| Applications | IF/IHC/IP/WB |
| Host / Source | Rabbit |
| Reactivity | Human/Mouse/Rat |

PRODUCT PROPERTIES

| | |
|---------------------|---|
| Clonality | Monoclonal |
| Conjugation | Unconjugated |
| Purification | Protein A/G purified from cell culture supernatant |
| Dilution Range | IF: 1:50-1:200 |
| Formulation | 0.01M PBS |
| Isotype | IgG |
| Storage Instruction | Suitable for storage at +4°C between 1-2 weeks. For longer term store at -20°C for up to 12 months. |

TARGET INFORMATION

| | |
|-------------|-----------------------------|
| Gene ID | 10018 |
| Gene Symbol | BCL2L11 |
| UniProt ID | B2L11_HUMAN |

ADDITIONAL INFORMATION

| | |
|------|--|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|--|