

**Anti-Recombinant-TCF21 antibody [RM1E27] (STJA0023021)****GENERAL INFORMATION**

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | FC/IP/WB           |
| Host / Source | Rabbit             |
| Reactivity    | Human              |

**PRODUCT PROPERTIES**

|                     |                                                                                                     |
|---------------------|-----------------------------------------------------------------------------------------------------|
| Clonality           | Monoclonal                                                                                          |
| Conjugation         | Unconjugated                                                                                        |
| Purification        | Protein A/G purified from cell culture supernatant                                                  |
| Dilution Range      | FCM: 1:20-1:100                                                                                     |
| Formulation         | 0.01M PBS                                                                                           |
| Isotype             | IgG                                                                                                 |
| Storage Instruction | Suitable for storage at +4°C between 1-2 weeks. For longer term store at -20°C for up to 12 months. |

**TARGET INFORMATION**

|             |                             |
|-------------|-----------------------------|
| Gene ID     | <a href="#">6943</a>        |
| Gene Symbol | <a href="#">TCF21</a>       |
| UniProt ID  | <a href="#">TCF21_HUMAN</a> |

**ADDITIONAL INFORMATION**

|      |                                                                                                                      |
|------|----------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|----------------------------------------------------------------------------------------------------------------------|