

**Anti-Recombinant-Phospho-RPS6KA5/MSK1-Ser376 antibody [RM1F58] (STJA0023151)****GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	IHC/IP/WB
Host / Source	Rabbit
Reactivity	Human/Rat

**PRODUCT PROPERTIES**

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein A/G purified from cell culture supernatant
Dilution Range	IHC: 1:50-1:100
Formulation	0.01M PBS
Isotype	IgG
Storage Instruction	Suitable for storage at +4°C between 1-2 weeks. For longer term store at -20°C for up to 12 months.

**TARGET INFORMATION**

Gene ID	<a href="#">9252</a>
Gene Symbol	<a href="#">RPS6KA5</a>
UniProt ID	<a href="#">KS6A5_HUMAN</a>

**ADDITIONAL INFORMATION**

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--