

Anti-Recombinant-AGR2/HPC8 antibody [RM1G83] (STJA0023275)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IF/IHC/IP/WB
Host / Source	Rabbit
Reactivity	Human/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein A/G purified from cell culture supernatant
Dilution Range	IF: 1:50-1:200
Formulation	0.01M PBS
Isotype	IgG
Storage Instruction	Suitable for storage at +4°C between 1-2 weeks. For longer term store at -20°C for up to 12 months.

TARGET INFORMATION

Gene ID	10551
Gene Symbol	AGR2
UniProt ID	AGR2_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--