

**Anti-Recombinant-Mono-methyl-Histone H2A-K118 antibody [RM1J80] (STJA0023470)****GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	IF/IHC/WB
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

**PRODUCT PROPERTIES**

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein A/G purified from cell culture supernatant
Dilution Range	IF: 1:50-1:200
Formulation	0.01M PBS
Isotype	IgG
Storage Instruction	Suitable for storage at +4°C between 1-2 weeks. For longer term store at -20°C for up to 12 months.

**TARGET INFORMATION**

Gene ID	<a href="#">3012/8335</a>
Gene Symbol	<a href="#">H2AC4.H2AC8</a>
UniProt ID	<a href="#">H2A1B_HUMAN</a>

**ADDITIONAL INFORMATION**

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--