

## Anti-Recombinant-TFAP2A antibody [RM1K36] (STJA0023525)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IP/WB
Host / Source	Rabbit
Reactivity	Rat

### PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein A/G purified from cell culture supernatant
Dilution Range	IP: 1:20
Formulation	0.01M PBS
Isotype	IgG
Storage Instruction	Suitable for storage at +4°C between 1-2 weeks. For longer term store at -20°C for up to 12 months.

### TARGET INFORMATION

Gene ID	<a href="#">7020</a>
Gene Symbol	<a href="#">TFAP2A</a>
UniProt ID	<a href="#">AP2A_HUMAN</a>

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--