

**Anti-Recombinant-Phospho-KRT8/CK-8-Ser23 antibody [RM1K43] (STJA0023532)****GENERAL INFORMATION**

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | IF/IHC/IP/WB       |
| Host / Source | Rabbit             |
| Reactivity    | Human              |

**PRODUCT PROPERTIES**

|                     |   |
|---------------------|---|
| Clonality           | Monoclonal  |
| Conjugation         | Unconjugated  |
| Purification        | Protein A/G purified from cell culture supernatant  |
| Dilution Range      | IF: 1:50-1:200  |
| Formulation         | 0.01M PBS   |
| Isotype             | IgG   |
| Storage Instruction | Suitable for storage at +4°C between 1-2 weeks. For longer term store at -20°C for up to 12 months. |

**TARGET INFORMATION**

|             |                            |
|-------------|----------------------------|
| Gene ID     | <a href="#">3856</a>       |
| Gene Symbol | <a href="#">KRT8</a>       |
| UniProt ID  | <a href="#">K2C8_HUMAN</a> |

**ADDITIONAL INFORMATION**

|      |  |
|------|--|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|--|