

Anti-Recombinant-PFKM antibody [RM1L51] (STJA0023639)

GENERAL INFORMATION

| | |
|---------------|--------------------|
| Product Type | Primary antibodies |
| Applications | FC/IF/WB |
| Host / Source | Rabbit |
| Reactivity | Human/Mouse/Rat |

PRODUCT PROPERTIES

| | |
|---------------------|---|
| Clonality | Monoclonal |
| Conjugation | Unconjugated |
| Purification | Protein A/G purified from cell culture supernatant |
| Dilution Range | FCM: 1:20-1:100 |
| Formulation | 0.01M PBS |
| Isotype | IgG |
| Storage Instruction | Suitable for storage at +4°C between 1-2 weeks. For longer term store at -20°C for up to 12 months. |

TARGET INFORMATION

| | |
|-------------|-----------------------------|
| Gene ID | 5213 |
| Gene Symbol | PFKM |
| UniProt ID | PFKAM_HUMAN |

ADDITIONAL INFORMATION

| | |
|------|--|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|--|