

**Anti-Recombinant-DLD antibody [RM1M21] (STJA0023708)****GENERAL INFORMATION**

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | IF/IHC/WB          |
| Host / Source | Rabbit             |
| Reactivity    | Human/Mouse/Rat    |

**PRODUCT PROPERTIES**

|                     |                                                                                                     |
|---------------------|-----------------------------------------------------------------------------------------------------|
| Clonality           | Monoclonal                                                                                          |
| Conjugation         | Unconjugated                                                                                        |
| Purification        | Protein A/G purified from cell culture supernatant                                                  |
| Dilution Range      | IF: 1:50-1:200                                                                                      |
| Formulation         | 0.01M PBS                                                                                           |
| Isotype             | IgG                                                                                                 |
| Storage Instruction | Suitable for storage at +4°C between 1-2 weeks. For longer term store at -20°C for up to 12 months. |

**TARGET INFORMATION**

|             |                            |
|-------------|----------------------------|
| Gene ID     | <a href="#">1738</a>       |
| Gene Symbol | <a href="#">DLD</a>        |
| UniProt ID  | <a href="#">DLDH_HUMAN</a> |

**ADDITIONAL INFORMATION**

|      |                                                                                                                      |
|------|----------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|----------------------------------------------------------------------------------------------------------------------|