

## Anti-Recombinant-CD324/CDH1 antibody [RM1N75] (STJA0023861)

### GENERAL INFORMATION

|               |                         |
|---------------|-------------------------|
| Product Type  | Primary antibodies      |
| Applications  | IF/IHC/IP/WB            |
| Host / Source | Rabbit                  |
| Reactivity    | Human/Mouse/Rat/Hamster |

### PRODUCT PROPERTIES

|                     |   |
|---------------------|---|
| Clonality           | Monoclonal  |
| Conjugation         | Unconjugated  |
| Purification        | Protein A/G purified from cell culture supernatant  |
| Dilution Range      | IF: 1:50-1:200  |
| Formulation         | 0.01M PBS   |
| Isotype             | IgG   |
| Storage Instruction | Suitable for storage at +4°C between 1-2 weeks. For longer term store at -20°C for up to 12 months. |

### TARGET INFORMATION

|             |                             |
|-------------|-----------------------------|
| Gene ID     | <a href="#">999</a>         |
| Gene Symbol | <a href="#">CDH1</a>        |
| UniProt ID  | <a href="#">CADH1_HUMAN</a> |

### ADDITIONAL INFORMATION

|      |  |
|------|--|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|--|