

Anti-Recombinant-Phospho-Src-Tyr530 antibody [RM1N77] (STJA0023863)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	IP/WB
Host / Source	Rabbit
Reactivity	Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein A/G purified from cell culture supernatant
Dilution Range	IP: 1:20
Formulation	0.01M PBS
Isotype	IgG
Storage Instruction	Suitable for storage at +4°C between 1-2 weeks. For longer term store at -20°C for up to 12 months.

TARGET INFORMATION

Gene ID	6714
Gene Symbol	SRC
UniProt ID	SRC_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--