

Anti-Recombinant-ACVR2A/ACTR-IIA antibody [RM1T20] (STJA0024301)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	WB
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein A/G purified from cell culture supernatant
Dilution Range	WB: 1:1000-1:2000
Formulation	0.01M PBS
Isotype	IgG
Storage Instruction	Suitable for storage at +4°C between 1-2 weeks. For longer term store at -20°C for up to 12 months.

TARGET INFORMATION

Gene ID	92
Gene Symbol	ACVR2A
UniProt ID	AVR2A_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--