

Anti-Recombinant-CD26/DPP4 antibody [RM1T33] (STJA0024314)

GENERAL INFORMATION		
	Product Type	Primary antibodies
	Applications	WB
	Host / Source	Rabbit
	Reactivity	Human/Mouse/Rat
PRODUCT PROPERTIES		
	Clonality	Monoclonal
	Conjugation	Unconjugated
	Purification	Protein A/G purified from cell culture supernatant
	Dilution Range	WB: 1:1000-1:2000
	Formulation	0.01M PBS
	Isotype	IgG
	Storage Instruction	Suitable for storage at +4°C between 1-2 weeks. For longer term store at -20°C for up to 12 months.
TARGET INFORMATION		
	Gene ID	1803
	Gene Symbol	DPP4
	UniProt ID	DPP4_HUMAN
ADDITIONAL INFORMATION		
	Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.