

**Anti-Recombinant-CASP1/Caspase-1 antibody [RM1T75] (STJA0024356)****GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	IP/WB
Host / Source	Rabbit
Reactivity	Human

**PRODUCT PROPERTIES**

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein A/G purified from cell culture supernatant
Dilution Range	IP: 1:20
Formulation	0.01M PBS
Isotype	IgG
Storage Instruction	Suitable for storage at +4°C between 1-2 weeks. For longer term store at -20°C for up to 12 months.

**TARGET INFORMATION**

Gene ID	<a href="#">834</a>
Gene Symbol	<a href="#">CASP1</a>
UniProt ID	<a href="#">CASP1_HUMAN</a>

**ADDITIONAL INFORMATION**

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--