

## Anti-Recombinant-ACTA1 antibody [RM2A85] (STJA0025059)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC/IF/WB
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat/Chicken

### PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein A/G purified from cell culture supernatant
Dilution Range	FCM: 1:20-1:100
Formulation	0.01M PBS
Isotype	IgG
Storage Instruction	Suitable for storage at +4°C between 1-2 weeks. For longer term store at -20°C for up to 12 months.

### TARGET INFORMATION

Gene ID	<a href="#">58</a>
Gene Symbol	<a href="#">ACTA1</a>
UniProt ID	<a href="#">ACTS_HUMAN</a>

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--