

Anti-Recombinant-Phospho-Histone H3-Ser10/Thr11 antibody [RM2B17] (STJA0025090)

GENERAL INFORMATION		
	Product Type	Primary antibodies
	Applications	IF/IHC/IP/WB
	Host / Source	Rabbit
	Reactivity	Human/Mouse/Rat
PRODUCT PROPERTIES		
	Clonality	Monoclonal
	Conjugation	Unconjugated
	Purification	Protein A/G purified from cell culture supernatant
	Dilution Range	IF: 1:50-1:200, IHC: 1:50-1:100, IP: 1:20, WB: 1:500-1:1000
	Formulation	0.01M PBS, pH 7.4, 0.05% BSA, 50% Glycerol, 0.05% Sodium Azide
	Isotype	IgG
	Storage Instruction	Suitable for storage at +4°C between 1-2 weeks. For longer term store at -20°C for up to 12 months.
TARGET INFORMATION		
	Gene ID	8350/8351/8352/8353/8354/8355/8356/8357/8358/8968
	Gene Symbol	H3C1.H3C2.H3C3.H3C4.H3C6.H3C7.H3C8.H3C10.H3C11.H3C12
	UniProt ID	H31_HUMAN
ADDITIONAL INFORMATION		
	Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.