

## Anti-Recombinant-NTHL1 antibody [RM2B35] (STJA0025108)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | IF/IHC/WB          |
| Host / Source | Rabbit             |
| Reactivity    | Human              |

### PRODUCT PROPERTIES

|                     |   |
|---------------------|---|
| Clonality           | Monoclonal  |
| Conjugation         | Unconjugated  |
| Purification        | Protein A/G purified from cell culture supernatant  |
| Dilution Range      | IF: 1:50-1:200, IHC: 1:100-1:200, WB: 1:1000-1:2000   |
| Formulation         | 0.01M PBS, pH 7.4, 0.05% BSA, 50% Glycerol, 0.05% Sodium Azide                                      |
| Isotype             | IgG   |
| Storage Instruction | Suitable for storage at +4°C between 1-2 weeks. For longer term store at -20°C for up to 12 months. |

### TARGET INFORMATION

|             |                           |
|-------------|---------------------------|
| Gene ID     | <a href="#">4913</a>      |
| Gene Symbol | <a href="#">NTHL1</a>     |
| UniProt ID  | <a href="#">NTH_HUMAN</a> |

### ADDITIONAL INFORMATION

|      |  |
|------|--|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|--|