

**Anti-Recombinant-Phospho-PRKCE-S729 antibody [RM2C13] (STJA0025185)****GENERAL INFORMATION**

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | WB                 |
| Host / Source | Rabbit             |
| Reactivity    | Human/Rat          |

**PRODUCT PROPERTIES**

|                     |   |
|---------------------|---|
| Clonality           | Monoclonal  |
| Conjugation         | Unconjugated  |
| Purification        | Protein A/G purified from cell culture supernatant  |
| Dilution Range      | WB: 1:1000-1:2000   |
| Formulation         | 0.01M PBS, pH 7.4, 0.05% BSA, 50% Glycerol, 0.05% Sodium Azide                                      |
| Isotype             | IgG   |
| Storage Instruction | Suitable for storage at +4°C between 1-2 weeks. For longer term store at -20°C for up to 12 months. |

**TARGET INFORMATION**

|         |  |
|---------|--|
| Gene ID | <a href="#">5581</a><br><a href="#">5579</a><br><a href="#">5578</a><br><a href="#">PRKCE</a><br><a href="#">PRKCB</a><br><a href="#">KPCE_HUMAN</a><br><a href="#">KPCB_HUMAN</a><br><a href="#">KPCA_HUMAN</a> |
|---------|--|

**ADDITIONAL INFORMATION**

|      |  |
|------|--|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|--|