

## Anti-Recombinant-UBE3A antibody [RM2C54] (STJA0025226)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | IF/IP/WB           |
| Host / Source | Rabbit             |
| Reactivity    | Human/Mouse/Rat    |

### PRODUCT PROPERTIES

|                     |                                                                                                     |
|---------------------|-----------------------------------------------------------------------------------------------------|
| Clonality           | Monoclonal                                                                                          |
| Conjugation         | Unconjugated                                                                                        |
| Purification        | Protein A/G purified from cell culture supernatant                                                  |
| Dilution Range      | IF: 1:50-1:200, IP: 1:20, WB: 1:500-1:1000                                                          |
| Formulation         | 0.01M PBS, pH 7.4, 0.05% BSA, 50% Glycerol, 0.05% Sodium Azide                                      |
| Isotype             | IgG                                                                                                 |
| Storage Instruction | Suitable for storage at +4°C between 1-2 weeks. For longer term store at -20°C for up to 12 months. |

### TARGET INFORMATION

|             |                             |
|-------------|-----------------------------|
| Gene ID     | <a href="#">7337</a>        |
| Gene Symbol | <a href="#">UBE3A</a>       |
| UniProt ID  | <a href="#">UBE3A_HUMAN</a> |

### ADDITIONAL INFORMATION

|      |                                                                                                                      |
|------|----------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|----------------------------------------------------------------------------------------------------------------------|