

## Anti-Recombinant-SRSF1 antibody [RM2D18] (STJA0025289)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IF/IHC/WB
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

### PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein A/G purified from cell culture supernatant
Dilution Range	IF: 1:50-1:200, IHC: 1:50-1:100, WB: 1:500-1:1000
Formulation	0.01M PBS, pH 7.4, 0.05% BSA, 50% Glycerol, 0.05% Sodium Azide
Isotype	IgG
Storage Instruction	Suitable for storage at +4°C between 1-2 weeks. For longer term store at -20°C for up to 12 months.

### TARGET INFORMATION

Gene ID	<a href="#">6426</a>
Gene Symbol	<a href="#">SRSF1</a>
UniProt ID	<a href="#">SRSF1_HUMAN</a>

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--