

Anti-Recombinant-CAMK1 antibody [RM2F15] (STJA0025484)

GENERAL INFORMATION

| | |
|---------------|--------------------|
| Product Type | Primary antibodies |
| Applications | WB |
| Host / Source | Rabbit |
| Reactivity | Mouse/Rat |

PRODUCT PROPERTIES

| | |
|---------------------|---|
| Clonality | Monoclonal |
| Conjugation | Unconjugated |
| Purification | Protein A/G purified from cell culture supernatant |
| Dilution Range | WB: 1:500-1:1000 |
| Formulation | 0.01M PBS, pH 7.4, 0.05% BSA, 50% Glycerol, 0.05% Sodium Azide |
| Isotype | IgG |
| Storage Instruction | Suitable for storage at +4°C between 1-2 weeks. For longer term store at -20°C for up to 12 months. |

TARGET INFORMATION

| | |
|-------------|-----------------------------|
| Gene ID | 8536 |
| Gene Symbol | CAMK1 |
| UniProt ID | KCC1A_HUMAN |

ADDITIONAL INFORMATION

| | |
|------|--|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|--|