

Anti-Recombinant-Phospho-CDC37-Ser13 antibody [RM2H03] (STJA0025670)

GENERAL INFORMATION		
	Product Type	Primary antibodies
	Applications	IP/WB
	Host / Source	Rabbit
	Reactivity	Human/Mouse/Rat
PRODUCT PROPERTIES		
	Clonality	Monoclonal
	Conjugation	Unconjugated
	Purification	Protein A/G purified from cell culture supernatant
	Dilution Range	IP: 1:20, WB: 1:500-1:1000
	Formulation	0.01M PBS, pH 7.4, 0.05% BSA, 50% Glycerol, 0.05% Sodium Azide
	Isotype	IgG
	Storage Instruction	Suitable for storage at +4°C between 1-2 weeks. For longer term store at -20°C for up to 12 months.
TARGET INFORMATION		
	Gene ID	11140
	Gene Symbol	CDC37
	UniProt ID	CDC37_HUMAN
ADDITIONAL INFORMATION		
	Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.