

Anti-Recombinant-LGR6 antibody [RM2M38] (STJA0026101)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IF/IHC/WB
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein A/G purified from cell culture supernatant
Dilution Range	IF: 1:50-1:200, IHC: 1:50-1:100, WB: 1:500-1:1000
Formulation	0.01M PBS, pH 7.4, 0.05% BSA, 50% Glycerol, 0.05% Sodium Azide
Isotype	IgG
Storage Instruction	Suitable for storage at +4°C between 1-2 weeks. For longer term store at -20°C for up to 12 months.

TARGET INFORMATION

Gene ID	59352
Gene Symbol	LGR6
UniProt ID	LGR6_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--