

## Anti-Recombinant-UBA2/SAE2 antibody [RM2N79] (STJA0026241)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | IF/IP/WB           |
| Host / Source | Rabbit             |
| Reactivity    | Human/Rat          |

### PRODUCT PROPERTIES

|                     |   |
|---------------------|---|
| Clonality           | Monoclonal  |
| Conjugation         | Unconjugated  |
| Purification        | Protein A/G purified from cell culture supernatant  |
| Dilution Range      | IF: 1:50-1:200, IP: 1:20, WB: 1:500-1:1000  |
| Formulation         | 0.01M PBS, pH 7.4, 0.05% BSA, 50% Glycerol, 0.05% Sodium Azide                                      |
| Isotype             | IgG   |
| Storage Instruction | Suitable for storage at +4°C between 1-2 weeks. For longer term store at -20°C for up to 12 months. |

### TARGET INFORMATION

|             |                            |
|-------------|----------------------------|
| Gene ID     | <a href="#">10054</a>      |
| Gene Symbol | <a href="#">UBA2</a>       |
| UniProt ID  | <a href="#">SAE2_HUMAN</a> |

### ADDITIONAL INFORMATION

|      |  |
|------|--|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|--|