

## Anti-Recombinant-CAMK2A antibody [RM2P05] (STJA0026266)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC/WB
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

### PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein A/G purified from cell culture supernatant
Dilution Range	FCM: 1:20-1:100, WB: 1:1000-1:2000
Formulation	0.01M PBS, pH 7.4, 0.05% BSA, 50% Glycerol, 0.05% Sodium Azide
Isotype	IgG
Storage Instruction	Suitable for storage at +4°C between 1-2 weeks. For longer term store at -20°C for up to 12 months.

### TARGET INFORMATION

Gene ID	<a href="#">816</a> <a href="#">818</a> <a href="#">815</a> <a href="#">CAMK2B</a> <a href="#">CAMK2G</a> <a href="#">KCC2B_HUMAN</a> <a href="#">KCC2G_HUMAN</a> <a href="#">KCC2A_HUMAN</a>
---------	--

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--