

## Anti-Recombinant-PPCS antibody [RM2RM90] (STJA0026549)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IF/IHC/WB
Host / Source	Rabbit
Reactivity	Human

### PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein A/G purified from cell culture supernatant
Dilution Range	IF: 1:50-1:200, IHC: 1:100-1:200, WB: 1:1000-1:2000
Formulation	0.01M PBS, pH 7.4, 0.05% BSA, 50% Glycerol, 0.05% Sodium Azide
Isotype	IgG
Storage Instruction	Suitable for storage at +4°C between 1-2 weeks. For longer term store at -20°C for up to 12 months.

### TARGET INFORMATION

Gene ID	<a href="#">79717</a>
Gene Symbol	<a href="#">PPCS</a>
UniProt ID	<a href="#">PPCS_HUMAN</a>

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--