

**Anti-Recombinant-Phospho-Erk1 + Erk2-T185/Y187 antibody [RM1B13] (STJA0026783)****GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	IF/IP/WB
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

**PRODUCT PROPERTIES**

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein A/G purified from cell culture supernatant
Dilution Range	IF: 1:50-1:200, IP: 1:20-1:50, WB: 1:1000-1:2000
Formulation	0.01M PBS, pH 7.4, 0.05% BSA, 50% Glycerol, 0.05% Sodium Azide
Isotype	IgG
Storage Instruction	Suitable for storage at +4°C between 1-2 weeks. For longer term store at -20°C for up to 12 months.

**TARGET INFORMATION**

Gene ID	<a href="#">5595</a> <a href="#">5594</a>
Gene Symbol	<a href="#">MAPK3</a> <a href="#">MAPK1</a>
UniProt ID	<a href="#">MK03_HUMAN</a> <a href="#">MK01_HUMAN</a>

**ADDITIONAL INFORMATION**

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--