

Anti-Recombinant-TUBB4B antibody [RM1B51] (STJA0026821)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC/IF/IHC/IP/WB
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein A/G purified from cell culture supernatant
Dilution Range	FCM: 1:20-1:100, IF: 1:50-1:200, IHC: 1:100-1:200, IP: 1:20-1:50, WB: 1:1000-1:2000
Formulation	0.01M PBS, pH 7.4, 0.05% BSA, 50% Glycerol, 0.05% Sodium Azide
Isotype	IgG
Storage Instruction	Suitable for storage at +4°C between 1-2 weeks. For longer term store at -20°C for up to 12 months.

TARGET INFORMATION

Gene ID [7280](#)
[10383](#)
[347733](#)
[TUBB2A](#)
[TUBB4B](#)
[TBB2A_HUMAN](#)
[TBB4B_HUMAN](#)
[TBB2B_HUMAN](#)

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.