

Anti-Recombinant-NRAS/Ras antibody [RM1B63] (STJA0026833)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	FC/IP/WB
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein A/G purified from cell culture supernatant
Dilution Range	FCM: 1:20-1:100, IP: 1:20-1:50, WB: 1:1000-1:5000
Formulation	0.01M PBS, pH 7.4, 0.05% BSA, 50% Glycerol, 0.05% Sodium Azide
Isotype	IgG
Storage Instruction	Suitable for storage at +4°C between 1-2 weeks. For longer term store at -20°C for up to 12 months.

TARGET INFORMATION

Gene ID	4893 3845 3265 NRAS KRAS RAS_N_HUMAN RASK_HUMAN RASH_HUMAN
---------	---

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--