

Anti-SLC1A5 antibody [SAA0181] (STJA0026967)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Human
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	0.58 mg/ml
Conjugation	Unconjugated
Purification	Protein A/G purified from cell culture supernatant.
Formulation	0.01M PBS, pH 7.4.
Isotype	IgG1k
Storage Instruction	Suitable for storage at +4°C between 1-2 weeks. For longer term store at -20°C for up to 12 months.

TARGET INFORMATION

Gene ID	6510
Gene Symbol	SLC1A5
UniProt ID	AAAT_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--