

Anti-IL1RL1/ST2 antibody [SAA0139] (STJA0027008)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC
Host / Source	Human
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/ml
Conjugation	Unconjugated
Purification	Protein A/G purified from cell culture supernatant.
Formulation	0.01M PBS, pH 7.4.
Isotype	IgG1
Storage Instruction	Suitable for storage at +4°C between 1-2 weeks. For longer term store at -20°C for up to 12 months.

TARGET INFORMATION

Gene ID	9173
Gene Symbol	IL1RL1
UniProt ID	ILRL1_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--