

## Anti-CD49c/ITGA3 antibody [84H7] (STJA0027311)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/FC/IP/WB
Host / Source	Mouse
Reactivity	Human

### PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/ml
Conjugation	Unconjugated
Purification	Protein A/G purified from cell culture supernatant.
Formulation	0.01M PBS, pH 7.4.
Isotype	IgG1k
Storage Instruction	Suitable for storage at +4°C between 1-2 weeks. For longer term store at -20°C for up to 12 months.

### TARGET INFORMATION

Gene ID	<a href="#">3675</a>
Gene Symbol	<a href="#">ITGA3</a>
UniProt ID	<a href="#">ITA3_HUMAN</a>

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--