

**Anti-Alpha 4 Beta 7/ITGA4 & ITGB7 antibody [SAA0020] (STJA0027312)****GENERAL INFORMATION**

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | FC                 |
| Host / Source | Human              |
| Reactivity    | Human              |

**PRODUCT PROPERTIES**

|                     |                                                                                                     |
|---------------------|-----------------------------------------------------------------------------------------------------|
| Clonality           | Monoclonal                                                                                          |
| Concentration       | 1 mg/ml                                                                                             |
| Conjugation         | Unconjugated                                                                                        |
| Purification        | Protein A/G purified from cell culture supernatant.                                                 |
| Formulation         | 0.01M PBS, pH 7.4.                                                                                  |
| Isotype             | IgG2k                                                                                               |
| Storage Instruction | Suitable for storage at +4°C between 1-2 weeks. For longer term store at -20°C for up to 12 months. |

**ADDITIONAL INFORMATION**

|      |                                                                                                                      |
|------|----------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|----------------------------------------------------------------------------------------------------------------------|