

Anti-CD8b/Lyt3 antibody [35.17.2#] (STJA0027419)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC/IF/IP/RIA/WB
Host / Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/ml
Conjugation	Unconjugated
Purification	Protein A/G purified from cell culture supernatant.
Formulation	0.01M PBS, pH 7.4.
Isotype	IgG1k
Storage Instruction	Suitable for storage at +4°C between 1-2 weeks. For longer term store at -20°C for up to 12 months.

TARGET INFORMATION

Gene ID	12526
Gene Symbol	Cd8b
UniProt ID	CD8B_MOUSE

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--