

Anti-CD11c/ITGAX antibody [N418] (STJA0027429)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	FC/IHC/IP/IF
Host / Source	Human
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	2.48 mg/ml
Conjugation	Unconjugated
Purification	Protein A/G purified from cell culture supernatant.
Formulation	0.01M PBS, pH 7.4.
Isotype	IgG1k
Storage Instruction	Suitable for storage at +4°C between 1-2 weeks. For longer term store at -20°C for up to 12 months.

TARGET INFORMATION

Gene ID	16411
Gene Symbol	Itgax
UniProt ID	ITAX_MOUSE

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--