

Anti-NOTCH1 antibody [SAA0270] (STJA0027502)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/FC
Host / Source	Rat
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/ml
Conjugation	Unconjugated
Purification	Protein A/G purified from cell culture supernatant.
Formulation	0.01M PBS, pH 7.4.
Isotype	IgG1k
Storage Instruction	Suitable for storage at +4°C between 1-2 weeks. For longer term store at -20°C for up to 12 months.

TARGET INFORMATION

Gene ID	18128
Gene Symbol	Notch1
UniProt ID	NOTC1_MOUSE

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	----------------------------------------------------------------------------------------------------------------------