

Anti-IL12A/IL-12 p35/NKSF1 antibody [SAA0380] (STJA0027559)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	ELISA/FC
Host / Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Concentration	1 mg/ml
Conjugation	Unconjugated
Purification	Protein A/G purified from cell culture supernatant.
Formulation	0.01M PBS, pH 7.4, 0.09% Sodium azide.
Isotype	IgG2ak
Storage Instruction	Suitable for storage at +4°C between 1-2 weeks. For longer term store at -20°C for up to 12 months.

TARGET INFORMATION

Gene ID	3592
Gene Symbol	IL12A
UniProt ID	IL12A_HUMAN

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	--