

Anti-IL13 antibody [SAA0383] (STJA0027562)

GENERAL INFORMATION		
	Product Type	Primary antibodies
	Applications	ELISA/FC
	Host / Source	Human
	Reactivity	Human
PRODUCT PROPERTIES		
	Clonality	Monoclonal
	Concentration	1 mg/ml
	Conjugation	Unconjugated
	Purification	Protein A/G purified from cell culture supernatant.
	Formulation	0.01M PBS, pH 7.4, 0.09% Sodium azide.
	Isotype	IgG4L
	Storage Instruction	Suitable for storage at +4°C between 1-2 weeks. For longer term store at -20°C for up to 12 months.
TARGET INFORMATION		
	Gene ID	3596
	Gene Symbol	IL13
	UniProt ID	IL13_HUMAN
ADDITIONAL INFORMATION		
	Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.