

Anti-CXCL10/IP-10 antibody [SAA0413] (STJA0027590)

GENERAL INFORMATION

| | |
|---------------|--------------------|
| Product Type | Primary antibodies |
| Applications | ELISA/FC |
| Host / Source | Rat |
| Reactivity | Human |

PRODUCT PROPERTIES

| | |
|---------------------|---|
| Clonality | Monoclonal |
| Concentration | 1 mg/ml |
| Conjugation | Unconjugated |
| Purification | Protein A/G purified from cell culture supernatant. |
| Formulation | 0.01M PBS, pH 7.4, 0.09% Sodium azide. |
| Isotype | IgG1k |
| Storage Instruction | Suitable for storage at +4°C between 1-2 weeks. For longer term store at -20°C for up to 12 months. |

TARGET INFORMATION

| | |
|-------------|-----------------------------|
| Gene ID | 3627 |
| Gene Symbol | CXCL10 |
| UniProt ID | CXL10_HUMAN |

ADDITIONAL INFORMATION

| | |
|------|--|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|--|