

## Anti-CD2/LFA-2/OX-34 antibody [OX-34] (STJA0027655)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | FC/IF/WB           |
| Host / Source | Mouse              |
| Reactivity    | Rat                |

### PRODUCT PROPERTIES

|                     |   |
|---------------------|---|
| Clonality           | Monoclonal  |
| Concentration       | 1 mg/ml   |
| Conjugation         | Unconjugated  |
| Purification        | Protein A/G purified from cell culture supernatant.   |
| Formulation         | 0.01M PBS, pH 7.4, 0.09% Sodium azide.  |
| Isotype             | IgG2ak  |
| Storage Instruction | Suitable for storage at +4°C between 1-2 weeks. For longer term store at -20°C for up to 12 months. |

### TARGET INFORMATION

|             |                         |
|-------------|-------------------------|
| Gene ID     | <a href="#">497761</a>  |
| Gene Symbol | <a href="#">Cd2</a>     |
| UniProt ID  | <a href="#">CD2_RAT</a> |

### ADDITIONAL INFORMATION

|      |  |
|------|--|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|--|